

Amjo Corp Background:

Amjo is a highly successful family business. **Now celebrating our tenth year**, we have become the World's Largest UV Phototherapy Distributor.

We're now in several distinct businesses. Our primary focus remains Home Phototherapy and we have become the **World's largest Home Phototherapy Distributor**.

Some of Amjo Corp's websites are:

www.HomePhotoTherapy.com
www.WoodsLight.com
www.ThermotexUSA.com
www.LightsForHealth.com
www.amjo.net
www.SADlight.com
www.BiomoveUSA.com
www.BeautySkinUSA.com
www.CatchaSmoker.com
www.Curatron.com

**USA TOLL FREE:
1-877-BUY-AMJO
1-877-289-2656**

Amjo Corp

PO Box 8304
West Chester, OH 45069
USA Toll Free: 1-877-BUY-AMJO
1-877-289-2656

Phone: 513-942-2770
Fax: 513-942-2771
info@amjo.net
www.HomePhotoTherapy.com

Home Phototherapy: Is it the right treatment for you?



Celebrating our Tenth Year
Your Home for Phototherapy
www.HomePhotoTherapy.com

**USA Toll Free:
1-877-BUY-AMJO**

Today It's UVB Narrow Band



The Panosol II UVB
Narrow Band Panel

UVB Narrow Band is used today to treat:

- Psoriasis
- Vitiligo
- Many forms of Eczema
- Mycosis Fungoides
- Many other skin challenges



You'll need to ask your doctor about phototherapy and if it is the right choice for you. For many with Vitiligo, Psoriasis or other skin challenges, sometimes the process of treatment is as troublesome as the disease.

Should you need in-office treatment, you, the patient may incur lost time from work or school, be burdened with traveling to and from a treatment facility and of course be faced with the expense and inconvenience of the treatment itself.

Now there's the Panosol II™ and the Foldalite III™ — Effective, Economical Home Phototherapy Systems.

These Ultraviolet systems bring the most advanced state of the art treatment into your own home. With your physician's prescription and guidance, you can now treat yourself in the privacy of your own home and on your own schedule. **What's more our UV Phototherapy Products are often insurance reimbursable.**



The Foldalite III requires minimal storage space, it folds to a mere 9" depth. Easy-open doors move and reposition effortlessly. Adjustable wheels for any floor surface.

www.HomePhotoTherapy.com

Amjo Corp is the World's Largest Home Phototherapy Distributor!

UV Home Phototherapy

Home Phototherapy Systems are here now! The tools your dermatologist uses are now available for you at home. Home systems such as the Panosol II and the Foldalite III use the same lamps and ballasts as those in professional systems.

Compare Home Phototherapy with the challenge of three drives a week to your dermatologist's office.

Imagine,

- Complete Privacy!
- No Travel Time!
- No Waiting!
- Ready when you are!



Panosol II, optional wings being opened for use.



National Biological Corporation Home UVB systems use high output ultraviolet-B (UVB) or UVB Narrow Band (UVBNB) lamps. These technologies have been used for many years to help control many skin disorders. UV remains a primary treatment by most dermatologists.

The lightweight Panosol II and its cousin the sixteen lamp Foldalite III run on standard house current and require no special circuits or wiring making them a breeze to use.

The CPT or Controlled Prescription Timer™ assures your physician control of your treatment regimen. The National Biological Corp's Controlled Prescription Timer offers unrivaled safety features including watchdog "fail safe" circuitry.

Q. What is Phototherapy?

A. Phototherapy, sometimes referred to as light therapy is the use of specific wavelengths of light that have been shown to be very effective. This light occurs naturally as a component of sunlight and is called ultraviolet light.

Q. What is UV Light?

A. The UV Light used in phototherapy is defined as UVB (Medium Wavelength) UV light or UVA (Long Wavelength UV Light) and today UVB Narrow Band is a very narrow band of wavelengths centered around 311 nm which seem to provide the optimal treatment with very little risk of sunburn.

Q. Is UV Phototherapy Safe?

A. At the time of writing there are no documented cases of any skin cancers attributable to UVB or UVB Narrow Band light. Excessive exposure may cause problems but even this is yet to be proven. Indications are that UVB and UVB Narrow Band are very safe when used within the guidelines provided by your doctor.

Q. Which form of UV should I use, UVA, UVB or UVB Narrow Band?

A. Probably UVB Narrowband for most skin challenges. UVA requires the use of Psoralen or Methoxalen taken internally or as a cream to make your skin sensitive to the UVA light. These drugs can introduce other problems such as nausea or other discomfort. PUVA (UVA and Psoralen) requires the patient wear UV blocking eyewear in the 24 to 48 hour period following treatment as the eyes also become extremely sensitive to UV light.

Q. How does phototherapy compare to other forms of treatment.

A. Typically other treatments require prescription steroids or systemic drugs. These alternatives are costly, probably not effective for long term use and have the potential for serious side effects.

Q. Is home phototherapy expensive?

A. Today many insurance companies will cover the majority of the expense. Amjo can help you throughout

the insurance claim process. We file insurance claims on your behalf and upon receipt of the payment from your insurance company and any co-pay amount from you we will ship a system immediately.

Q. What's the warranty on Home Phototherapy Systems from Amjo?

A. The Panosol II Panels, our Hand/Foot II and of course the Foldalite III all come with a limited Three Year in-home warranty, excluding lamps which are warranted for ninety days. Our Handheld wands come with a one year limited warranty.

Selected Pricing* Information:

Foldalite III Full Body System

FLUNB16 - 16 Lamp UVB Narrow Band	\$6850.00 **
FLUNB32 - 32 Lamp UVB Narrow Band	Please Call

Six Foot Panels:

P2UNB64 - UVB Narrow Band (4 Lamp)	\$1990.00 **
P2UNB66 - UVB Narrow Band (6 Lamp)	\$2390.00 **
P2UNB66RW - Above with reflective wings	\$2790.00 **
P2UNB3D - UVB Narrow Band 10 Lamp	\$3990.00 **
P2UVB64 - UVB Broadband (4 Lamp)	\$1975.00 **
P2UVB64RW - Above with reflective wings	\$2395.00 **
P2UVA66 - UVA (6 Lamp)	\$1975.00 **

Two Foot Tall Panels:

P2UNB24 - UVB Narrow Band (4 Lamp)	\$1595.00 **
P2UNB28 - UVB Narrow Band (8 Lamp)	\$1890.00 **
P2UVB24 - UVB Broadband (4 Lamp)	\$1590.00 **

Hand/Foot II Units:

HFUNB28 - UVB Narrow Band	\$1990.00 **
HFUVB28 - UVB Broadband	\$1890.00 **
HFUVA28 - UVA Unit	\$1890.00 **

Handisol Series Handheld Wands

HSUNB - UVB Narrow Band Wand	\$595.00 **
HSUVB - UVB Broadband Wand	\$690.00 **
HSUVA - UVA Handheld Wand	\$690.00 **

* As always, pricing is subject to change.

** Shipping costs included in USA lower 48 states



The Panosol II 3D at the left will feature a total of ten lamps.

Will be priced at \$3,990.00 **
Delivery Incl.
 (USA lower 48 States)
 (Available in Late 2008)

Without a doubt the world's most popular hand held phototherapy wand is the Handisol NB or Dermalight 80 pictured below. Thousands of this "magic" wand have been shipped worldwide. The



Handisol NB UVB Narrow Band Phototherapy Wand.

psoriasis. For all skin challenges, the comb can be used as a spacer to fix the distance to your skin for consistent day to day treatment.

Whatever your choice, for Home Phototherapy, Be Sure to Call Amjo Today!

Physician Approval:

I authorize patient _____ to purchase a [UVA] [UVB] [UVB Narrow Band] (Please circle one) Phototherapy Device. This patient has been instructed to consult me on a regular basis for follow-up exams and is aware of the treatment procedures for this product.

Signature: _____

Physician Information:

Name: _____ Phone: _____

Street: _____

City, State, Zip: _____

Caution: USA Federal law restricts this device to sale by or on the order of Physician (USA Only)